



#### Special characteristics

- Easy and safe handling by a single adjusting knob
- Modern controller with clear-text display
- Electronic current limiting and safety switching, which avoids electric shocks
- Extensive aluminium plate electrodes for a homogeneous current density and to avoid allergies
- Mains operation, as battery operation would be a lot more expensive
- Stable suitcase, which also serves as basins for the hand or feet treatment

## idromed® 5 GS / PS

The effective iontophoresis appliances for the treatment of hyperhidrosis

#### Direct current (GS) or pulsed direct current (PS)

The direct current and pulsed direct current iontophoresis-therapy is an internationally established and scientifically acknowledged treatment procedure, which is known for years.

In scientific studies the therapy with pulsed direct current with the idromed® 5 PS is rated as the optimum therapy method, thanks to its new current type, especially for very sensitive patients, children and young people. The therapy current is almost not perceived. Possible side effects are clearly reduced. The hands or feet can be taken out of the basins during the therapy without the usually known "electric fence effect" (electric shock).



**Optional sponge electrodes for armpits** Art.-no. 995100



Art.-no. idromed® 5 GS: 100 000  
Art.-no. idromed® 5 PS: 110 000

#### Technical Data

Dimensions (L x W x H):  
290 x 390 x 110 mm  
Dimensions of the unit  
(L x W x H): 210 x 120 x 70 mm  
Total weight: approx. 2,5 kg  
(unit, suitcase and accessories)  
Power supply:  
230V / 50Hz / 20VA  
Subject to technical alterations